



VPro+TM

QUICK START GUIDE




YOUR BEST ORTHODONTIC EXPERIENCE YET

The VPro+ is highly portable and designed with your busy lifestyle in mind. Use for just 5 minutes each day, or as otherwise directed by your doctor, to compliment your orthodontic treatment.

SELFIES SPREAD SMILES
#propelortho


WHAT'S IN THE BOX?

+MOUTHPIECE




+OSCILLATOR

+CHARGING CASE



+MICRO-USB CABLE



+WALL ADAPTER



WHERE TO START?



CHARGE

Be sure to charge your VPro+ fully before first use.



DOWNLOAD

Download the VPro+ Mobile App to start tracking your usage.

App use is optional



CONNECT

Secure the Mouthpiece to the Oscillator (arrow points toward button).



BEGIN

Begin your first 5 minute session. Wear your orthodontics during use unless otherwise directed.



ALL CHARGED UP

To charge the device, remove the Mouthpiece and place the Oscillator in the Charging Case while it's connected to a power source. A green light indicator will illuminate to confirm a successful connection. Do not charge for more than 4 hours at once. Recommended charging time is once per week or as needed.



PAIRING IS CARING

Pair your new VPro+ with your phone to track device usage, share results with your treating clinician, and locate the device as needed. Search "VPro+" in the Apple App store or Google Play store and download. Follow the on-screen onboarding instructions. App use is optional for device operation.



CLEAN AS A WHISTLE

Your Mouthpiece should be cleaned after each use. Detach the Mouthpiece from the Oscillator. Rinse the Mouthpiece under warm water and let dry. Clean the Oscillator with warm water as needed and let dry prior to storage/charging.



HANDLE WITH CARE

The device should be stored in the charging case. Proper care should be taken to ensure that the device is not damaged. Devices should be stored at room temperature, avoid any extreme heat.

ALL OF THE LIGHTS



The VPro+ will display light notifications from the device button, unless otherwise noted.

Notification	Meaning
Blue (3 Blinks)	Session Start
Blue (Slow Blinking)	Charging in Progress
Green LED (1 Blink)	App has Connected
Green (3 Blinks)	Session End
Green (Slow Blinking)	Charging Complete
Green (Blinks and Vibrates)	Device Locate
Yellow LED (2 Blinks)	Successful App Pairing
Red (1 Blink)	Critical Battery
Red (2 Blinks)	Device Reset/Reboot
Red (Back of Device)	Firmware Update/Charging

FAQs

Do I need to use the VPro+ app?

No. You can successfully utilize your VPro+ device without downloading the mobile app. The VPro+ app conveniently allows you to track your device usage and share that information with your doctor or others.

Can I use the VPro+ more than once a day?

While the suggested use is one 5 minute session per day, there are no contra-indications to utilizing the device multiple times a day. Please consult your doctor.

How much pressure should I apply?

Bite down gently on the mouthpiece during use. There is no need to bite down harder than as needed to support the device.

I got a new VPro+ because my old one was lost/stolen. How do I pair my new device?

First, unpair your old device by navigating to settings in the VPro+ App. Select the “Reset Device” button. Confirm your selection and your device is now unpaired. Now you may follow the in-app onboarding screens.

What should I do if the VPro+ will not charge?

If the VPro+ LED indicator does not light up red when placed in the case, please check the case’s power connection. If the case is properly plugged into a power source using the provided VPro+ USB cord, please contact Propel Customer Service for further assistance.

QUESTIONS? COMMENTS?

You Talk. We Listen.
Contact us today with your questions,
comments, and feedback.

Call Toll-Free
(855) 377-6735

info@propelortho.com

propelortho.com/contact

